

## 1.0 Example: Add the Material Aluminum from the Library.




Figure 1-1 Adding the Material Aluminum from the Library.

Start the **Visual Editor**.

Click on **Data...Materials** to open the Material Panel.

Click on **Library** on the Material Panel menu.

In the Material Library Panel window, **click** on **aluminum**.

On the Material Library Panel menu, **click** on **Add**.

In the Material Panel, **verify** that **aluminum now appears** as a material that may be selected in the material list at the top of the panel.

On the Visual Editor Main Menu, **click** on **File...Save** and enter **ilibmat** as the file name.

